**Altura™ Endograft System**  
Instructions for Use

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## ENGLISH

READ ALL INSTRUCTIONS CAREFULLY. FAILURE TO PROPERLY FOLLOW THE INSTRUCTIONS, WARNINGS, AND PRECAUTIONS MAY LEAD TO SERIOUS SURGICAL CONSEQUENCES OR INJURY TO THE PATIENT.

### 1.0 Product Description

The Altura™ Endograft System is comprised of two (2) bilateral aortic endografts and two (2) iliac endografts. Each endograft is pre-loaded into a 14 F delivery catheter. Each catheter is intended to be introduced into the femoral artery over a guidewire under fluoroscopic guidance. The two aortic endografts are advanced into the abdominal aorta, aligned, and positioned. The endografts are deployed such that the proximal edge of the graft material is below the renal arteries within the non-aneurysmal neck. The aortic delivery catheters are removed over the indwelling guidewires. Two delivery catheters containing iliac endografts are then delivered over the guidewires and positioned with respect to the intended sealing zone. After the iliac endografts are deployed the delivery catheters are removed. Balloon catheters are then inserted over the guidewires to dilate the endografts. A completion aortogram is performed to verify positioning and assess for endoleaks. **1.1. Device Components**

#### 1.1.1. Endograft

The aortic (**Figure 1** and **Figure 2**) and iliac endografts (**Figure 3**) are fabricated from a braided nitinol (nickel titanium alloy) wire stent and an ultrathin, high density, polyester graft, which is affixed to the stent with sutures at the proximal and distal ends of the graft. The proximal segment of the aortic endograft is not covered by graft material and is for suprarenal fixation. The two D-shaped aortic endografts are positioned side by side so that the flat surfaces appose each other to form a circular graft cross-section. The two aortic endografts can be longitudinally aligned along the flat surface with an offset of up to 10 mm to accommodate the different anatomic positions of the left and right renal arteries. An iliac endograft can be inserted within each aortic endograft to form a complete endograft construct.

**Note:** The diameters and lengths discussed below and shown in **Table 1** and **Table 2** are approximate and for unconstrained endografts. Endografts will have different dimensions dependent on the native anatomy.

The unconstrained polyester graft of the aortic endograft is 90 mm long and the uncovered suprarenal stent is 25 mm long. The aortic endograft is designed to terminate in the distal aorta. The aortic endograft has three radiopaque markers on the flat surface that are used to determine if the aortic endografts properly appose each other. These three radiopaque markers are referred to as the septal radiopaque markers. The aortic endograft also has a renal radiopaque marker immediately proximal to the leading edge of the graft material. The suprarenal stent portion of the aortic endograft flares to accommodate the suprarenal aorta. Two anchors are positioned on the suprarenal stent of each aortic endograft to secure the device within the aorta. The diameter of the leg of the aortic endograft is 13 mm.

The aortic endograft has three configurations, as shown below in **Table 1**.

**Table 1. Aortic Endograft Configurations**

Catalog Number	Aortic Endograft Diameter (mm)	Non-aneurysmal Neck Treatment Range (mm)	Suprarenal Stent Diameter (mm)	Aortic Leg Diameter (mm)
01-AA-24-090	24	18-22	29	13
01-AA-27-090	27	21-25	33	
01-AA-30-090	30	24-28	36	

The unconstrained polyester graft of the iliac endograft is 65 mm long. The uncovered distal stent is 2 mm long. The iliac endograft has two radiopaque markers on the proximal end of the stent that are used to verify that the endograft has released from the delivery catheter after deployment. The iliac endograft is available in three sizes: 13, 17, and 21 mm. These sizes refer to the diameter of the trailing end of the endograft. For the 17 mm and 21 mm size iliac endografts, the iliac endograft diameter section is 20 mm long and is referred to as the iliac sealing zone. All three sizes of iliac endograft have an iliac leg diameter of 13 mm. The iliac leg diameter is the leading end of the endograft that connects to and seals with the aortic endograft.

The iliac endograft has three configurations, as shown below in **Table 2**.

**Table 2. Iliac Endograft Configurations**

Catalog Number	Iliac Endograft Diameter (mm)	Iliac Treatment Range (mm)	Iliac Leg Diameter (mm)
01-AI-13-065	13	8-11	13
01-AI-17-065	17	11-15	
01-AI-21-065	21	15-18	

The Altura Endograft System is constructed of an iliac and aortic delivery system and endograft prosthesis. The endograft prostheses are constructed of the following materials.

**Table 3. Endograft Materials**

Component	Material
Graft	Polyester
Stent, Anchors, Sleeve	Nitinol
Radio Opaque Markers	Tantalum
Suture	Ultra High Molecular Weight Polyethylene (UHMWPE)

#### 1.1.2. Delivery Catheter

The Altura Endograft System delivery catheters are single use disposable devices used to deliver the aortic (Figure 4 and Figure 5) and iliac (Figure 6 and Figure 7) endografts. The aortic endografts are deployed in a proximal to distal fashion while the iliac endografts are deployed in a distal to proximal fashion.

Both types of endografts are delivered using 14 F introducer sheath compatible delivery catheters that contain 0.035" (0.89 mm) guidewire compatible lumens and radiopaque atraumatic tips. Endograft deployment is accomplished by rotating the rear grip of the delivery catheter handle as indicated on the handle.

Following initial exposure, the aortic endograft is expanded by pulling back the stent expander. When the endograft is properly positioned and expanded, the suprarenal stent is released by removing the safety collar and advancing the stent release forward. Deployment of the aortic endograft is completed by rotating the rear grip. The iliac endograft is delivered by rotating the rear grip of the delivery catheter until the endograft is released.

### Intended Purpose

The Altura Endograft System is intended to exclude aneurysms from blood circulation in patients diagnosed with abdominal aortic aneurysm (AAA). It is intended for use only by suitably trained physicians who are experienced in the diagnosis and endovascular treatment of aneurysmal disease. Standard techniques for the use of vascular access sheaths, angiography, guidewires and contrast media should be employed.

### 2.0 Indications for Use

The Altura Endograft System is indicated for treatment of patients with abdominal aortic and/or aorto-iliac aneurysms having vascular morphology suitable for endovascular repair, including:

- Adequate iliac/femoral access that is compatible with vascular access techniques, devices, and/or accessories
- Iliac/femoral access morphology compatible with 14 F introducer sheaths
- Non-aneurysmal proximal aortic neck between the renal arteries and the aneurysm with:
  - Length of  $\geq 15$  mm
  - Inner wall diameter of 18-28 mm
  - Neck angulation of  $\leq 60^\circ$
- Distal common iliac landing zone with:
  - Inner wall diameter of 8-18 mm
  - Length of  $\geq 15$  mm

Patients participating in a clinical research study using the Altura Endograft System should refer to the specific Investigational Plan for study specific indications for use and eligibility criteria.

### 3.0 Contraindications

The Altura Endograft System is contraindicated in the following patients.

- Patients with a systemic infection who may be at increased risk of endovascular graft infection.
- Patients with known sensitivities or allergies to polyester, nitinol, tantalum, or polyethylene.

### 4.0 Warnings and Precautions

#### 4.1. General

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.
- The Altura Endograft System should only be used by physicians and teams trained in endovascular interventional and surgical techniques. Specific training expectations are described in **Section 9.1, Physician Training and Experience**.
- Always have a qualified vascular surgery team available during the implantation procedure in the event that conversion to an open surgical repair is necessary.

- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing clinical sequelae, an enlarging aneurysm, and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.

#### 4.2. Patient Selection, Treatment, and Follow-Up

- The Altura Endograft System is designed to treat aortic neck diameters between 18 mm and 28 mm, inclusive. The Altura Endograft System is designed to treat proximal aortic necks (distal to the lowest renal artery) of at least 15 mm in length. Iliac artery distal fixation site greater than or equal to 15 mm in length and diameters between 8 mm and 18 mm is required. These sizing measurements are critical to the performance of the endovascular repair.
- Key anatomical elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (>60 degrees) and irregular or circumferential thrombus and/or calcification at the sealing zones that may compromise sealing and fixation. In the presence of anatomical limitations, a longer seal zone may be required to obtain adequate sealing and fixation.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measure inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease, and/or calcification) should be compatible with vascular access techniques and 14 F vascular introducer sheaths. Vessels that are significantly calcified, occlusive, tortuous, or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization. A vascular conduit technique may be necessary to achieve success in some patients.
- The Altura Endograft System is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endograft.
- The Altura Endograft System is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvis/bowel ischemia.
- Multiple large, patent lumbar arteries, and a patent inferior mesenteric artery may predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The Altura Endograft System has not been evaluated in the following patient populations:
  - Traumatic aortic injury
  - Leaking, dissecting, pending rupture, or ruptured aneurysms
  - Mycotic aneurysms
  - Pseudoaneurysms resulting from previous graft placement
  - Revision of previously placed endovascular grafts or previous AAA surgery
  - Uncorrectable coagulopathy
  - Dialysis or compromised renal and/or hepatic function
  - Genetic connective tissue disease (e.g. Marfans or Ehlers-Danlos Syndromes)
  - Concomitant thoracic aortic or thoracoabdominal aneurysms
  - Active systemic infections
  - Pregnant or nursing females
  - Morbidly obese patients
  - Less than 18 years of age
  - Patients with less than 15 mm in length or greater than 60 degrees angulation of the proximal aortic neck.
- Successful patient selection requires specific imaging and accurate measurements; please see Section 4.3, Pre-Procedure Measurement Techniques and Imaging.

#### 4.3. Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging of a maximum reconstruction thickness of 3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenosis.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Altura Endograft System. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning the x-ray C-arm during angiography such that the origins of the renal arteries, and particularly the lowest patient renal artery, are well demonstrated prior to deployment of the proximal edge of the endograft material. Additionally, angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are well defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg components.

## Diameters

Utilizing CTA, diameter measurements should be determined by the vessel lumen to help with proper device sizing and device selection. Adjustment to the measured lumen diameters may be necessary in the presence of thrombus, plaque, or calcium that may remodel after implantation of the endografts. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral bifurcations at an axial thickness slice of 2 mm or less.

## Lengths

Utilizing CTA, length measurements should be determined to accurately assess the following: proximal infrarenal neck length, renal to aortic bifurcation length, common iliac length, iliac fixation length, and total length from renal to internal iliac artery length.

### 4.4. Device Selection

Strict adherence to the Altura Endograft System IFU sizing guide is strongly recommended when selecting the appropriate device size (**Table 1** and **Table 2**). Appropriate device oversizing has been incorporated in the IFU sizing guide. Sizing outside of the range can result in endoleak, fracture, migration, or lumen loss.

### 4.5. Implant Procedure

- Appropriate procedural imaging is required to successfully position the Altura Endograft System and assure accurate apposition to the vessel walls.
- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endograft and delivery catheter during preparation and insertion to decrease the risk of endograft contamination and infection.
- Do not bend or kink the delivery catheters prior to insertion. Doing so may damage the delivery catheters and endografts.
- Maintain guidewire position during delivery catheter insertion.
- Fluoroscopy should be used during introduction, tracking, positioning, deployment, repositioning, and removal to confirm proper operation of the delivery catheter components, proper placement of the endograft, and desired procedural outcome. The use of the Altura Endograft System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise (e.g. adequate hydration).
- Do not fully implant a single aortic endograft within the aorta without a paired aortic endograft positioned, aligned, and partially expanded alongside. Fully implanting only one aortic endograft at a time can result in misalignment and improper positioning of the endograft within the aortic neck, resulting in an inadequate seal.
- As the sheath and/or guidewire is withdrawn, anatomy and endograft position may change. Constantly monitor endograft position and perform angiography to check position as necessary.
- To avoid any twist in the endovascular endograft, during any rotation of the delivery catheter ensure rotation of the delivery catheter handle is translated into rotation of the distal end of the delivery catheter. Relieve built-up torque in the delivery catheter by twisting the delivery catheter handle in the opposite direction after rotating the distal end into the desired orientation.
- The Altura Endograft System incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional and angiographic devices in the region of the suprarenal stent.
- Do not rotate the aortic delivery catheter once the suprarenal stent is released.
- Proper alignment of the flat surface of the aortic endografts is critical to ensure a good seal across the septum of the endografts and the neck of the aorta. The flat surfaces of the aortic endografts must appose each other. The aortic endografts should not be offset by more than 10 mm.
- Unless medically indicated, do not deploy the Altura Endograft System in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Vessel occlusion may occur.
- Do not attempt to re-sheath the endograft by rotating the rear grip of the delivery catheter handle in the opposite direction as indicated after partial or completed deployment.
- Repositioning the endograft distally after deployment of the suprarenal stent may result in damage to the endograft, reduction in endograft lumen, and/or vessel injury.
- Inaccurate placement and/or incomplete sealing of the Altura Endograft System within the vessel may result in increased risk of endoleak, migration, or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Altura Endograft System may result in increased risk of migration of the endograft. Incorrect deployment or migration of the endograft may require surgical intervention.
- Verify iliac endograft size prior to insertion into patient. Once the suprarenal stent is released within the aorta, care must be taken to not apply tension to the aortic endograft as it is released from the delivery catheter. This can result in endograft elongation, lumen reduction, and misalignment.
- During removal, care must be taken to prevent catching the delivery catheter on the endograft, anatomy, or the vascular introducer sheath. Using a vascular introducer sheath other than the brand and size specified in **Section 9.3, Required Materials**

may result in difficulty in removing the delivery catheters or vessel damage. Do not use excessive force to advance or withdraw the delivery catheters when resistance is encountered.

- Care must be taken when advancing the iliac endograft through the previously deployed aortic endograft to prevent the iliac delivery catheter from catching on the aortic endograft.
- The iliac endograft should overlap at least 20 mm within the aortic endograft to ensure a seal between the iliac and aortic endografts. If a minimum 20 mm overlap is not achieved, a 13 mm iliac endograft should be placed over the junction.
- Care must be taken when deploying the iliac endograft to avoid unplanned internal iliac artery occlusion.
- Do not reposition the iliac endograft distally once the endograft makes contact with the iliac artery.
- Do not rotate the iliac delivery catheter once the iliac endograft has engaged the anatomy.
- Ensure that the top cap radiopaque marker of the iliac delivery catheter is outside the vascular introducer sheath before the iliac endograft is deployed.
- When injecting contrast into the iliac or aortic delivery catheters, injector pressure should not exceed 300 psi.
- The tip of the guidewire should remain in position within the thoracic aorta throughout the procedure.
- Use caution during manipulation of catheters, wires, endografts, and sheaths within the vasculature. Significant disturbances may dislodge fragments of plaque or thrombus, which can cause distal embolization or rupture of the aneurysm.
- Inadvertent partial deployment, misalignment, or migration of the endograft may require surgical intervention.
- Avoid damaging the endograft or disturbing endograft positioning after placement in the event reinstrumentation (secondary intervention) of the endograft is necessary.
- Ensure introducer sheath position is maintained during the entire implant procedure.
- Do not continue advancing any portion of the delivery catheter if resistance is felt during advancement of the guidewire or delivery catheter. Stop and assess the cause of resistance using fluoroscopic guidance; vessel, catheter, or endograft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis, or in calcified or tortuous vessels.
- Do not use a delivery catheter if the packaging is damaged, the sterile barrier is broken prior to opening the packaging, or the delivery catheter appears damaged after removal from the packaging.

#### 4.6. Compliant Balloon Use

- Do not inflate the balloon in the vessel outside of the endograft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care in inflating the balloon within the endograft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.

#### 4.7. MRI Information

Non-clinical testing has demonstrated that the Altura Endograft System endografts are MR Conditional. Patients with this device can be scanned safely in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0-Tesla (3.0T)
- Maximum spatial gradient field less than or equal to 30 T/m (3000 gauss/cm)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the Altura Endograft System endografts are expected to produce a maximum temperature rise of 6.47°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 11 mm for the Altura Endograft System endografts when imaged with a gradient echo or spin echo pulse sequence in a 3.0 Tesla MRI system.

## 5.0 Potential Adverse Events

Adverse events that may occur and/or require intervention include, but are not limited to:

- Allergic reaction and/or intolerance to contrast or procedural/post-procedural medications
- Amputation
- Anesthesia complications and subsequent attendant problems (e.g. aspiration)
- Aneurysm enlargement
- Aneurysm rupture
- Aortic damage and subsequent attendant problems (e.g., bleeding, dissection, perforations, rupture)
- Aorto-enteric fistula
- Arterial and venous thrombosis and/or pseudoaneurysm
- Arterial venous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications and subsequent attendant problems (e.g., ileus, transient ischemia, infarction, necrosis)

- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock and lower limb)
- Conversion to open surgery
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endograft: improper endograft placement; incomplete endograft deployment; endograft migration; separation of graft material from stent; occlusion; allergic reaction to endograft material, infection; stent fracture; graft material wear; graft material failure; dilatation; erosion; puncture of graft material and perigraft flow; endograft deformity and/or kinking, corrosion
- Endoleak
- Fever or localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection, impotence)
- Impotence
- Incomplete aneurysm repair
- Increased risk of cancer due to radiation exposure
- Infection of the aneurysm, access site, including abscess formation, transient fever and pain
- Life-time annual requirement for follow-up examinations (involving radiation exposure and intravenous contrast agents) to verify that the endograft is continuing to function
- Liver failure
- Lymph fistula / complications
- Mechanical failure of the device
- Neurologic damage, local or systemic and subsequent attendant problems (e.g., stroke, paraplegia, paraparesis, transient ischemia, paralysis)
- Occlusion of the endograft or native vessels
- Pain
- Pulmonary or respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Reduced blood supply to distal limb and pelvic organs and subsequent attendant problems (e.g., leg pain, edema, amputation, buttock pain, reduced sexual function)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, hypertension)
- Surgery
- Vascular spasm or vascular trauma and subsequent attendant problems (e.g., ilio-femoral vessel dissection, infection, bleeding, rupture, pain, hematoma, arteriovenous fistula)
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., infection, dehiscence, bleeding, pseudoaneurysm)

## 6.0 Patient Selection and Treatment

### Individualization and Treatment

Lombard Medical recommends that the Altura Endograft System endograft diameters be selected as describe in Table 1 and Table 2. All device diameters potentially necessary to complete the procedure should be available. Considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g. cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the Altura Endograft System
- Ability to tolerate general, regional, or local anesthesia
- Iliofemoral access vessel size and morphology (minimal thrombus, calcium, and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 14 French vascular introducer sheath
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
  - with a length of at least 15 mm,
  - with a lumen diameter no greater than 28 mm and no less than 18 mm,

- with an angle less than or equal to 60 degrees, and Iliac artery distal fixation site greater than or equal to 15 mm in length and lumen diameters 8 to 18 mm
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft the final treatment decision is at the discretion of the physician and patient.

## 7.0 Patient Counseling

The physician and patient (and/or family members) should review the risks and benefits when discussing the endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- The long-term performance of the Altura endografts has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g. endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements are required and should be considered a lifelong commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise each patient that it is important to seek prompt medical attention if they experience signs of limb occlusion, aneurysm enlargement, or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal, or groin pain; dizziness; fainting; rapid heartbeat, or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant. Men who undergo endovascular or open surgical repair may experience impotence.

Procedure-related risks include cardiac, pulmonary, neurologic, bowel, and bleeding complications. Device-related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture, and death. **See Section 5, Potential Adverse Events.**

## 8.0 How Supplied

- The Altura Endograft System is sterilized by ethylene oxide gas and provided in a peel-open package.
- The Altura Endograft System endograft is preloaded in the delivery catheter.
- The device is intended for single use only. Do not re-sterilize the device.
- The product is sterile if the package is unopened or undamaged. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken.
- Prior to use, verify the correct devices (quantity and size) have been supplied for the patient.
- Do not use after the "USE BY" (expiration) date printed on the label.
- Store in a cool, dry place.

## 9.0 Clinical Use Information

### 9.1. Physician Training and Experience

**CAUTION: Always have a qualified surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.**

**CAUTION: The Altura Endograft System should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Altura Endograft System are outlined below:**

Patient Selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair

- Knowledge of radiographic image interpretation, device selection, planning, and sizing A multidisciplinary team that has combined procedural experience with:
- Femoral cutdown, arteriotomy, and repair
- Percutaneous access and closure techniques
- Non-selective and selective guidewire and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

### 9.2. Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. Prior to use, verify correct devices (quantity and size) have been supplied for the patient.

### 9.3. Required Materials

- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Various syringes for flushing and balloon inflation (e.g., 10, 20, and 30 ml (cc))
- Heparinized saline solution or equivalent
- Sterile gauze pads
- Vascular access kit
- COOK® Medical 14.0 F introducers
- Assorted 0.035" (0.89 mm) guidewires of adequate length
- Stiff 0.035" (0.89 mm) guidewires, 260 or 300 cm length
- Angiographic flush catheter with radiopaque markers
- Compliant vascular occlusion balloons
- Power injector
- Various vascular introducers

### 9.4. Recommended Materials

- Various diameter dilators
- Vascular closure devices
- Snares

### 9.5. Device Sizing Guidelines

The choice of diameter should be determined from **Table 1** and **Table 2**. Undersizing or oversizing may result in incomplete sealing or compromised flow. Altura recommends that the Altura Endograft System endograft component diameters be selected as described in **Table 1** and **Table 2**. All device diameters potentially necessary to complete the procedure should be available.

#### Summary of Safety and Clinical Performance (SSCP)

Information produced by Lombard Medical under the Medical Device Regulation (EU) 2017/745. The Summary of Safety and Clinical Performance is intended to provide public access to an updated summary of clinical data and other information about the safety and clinical performance of the medical device.

SSCP information is located on the company website [www.lombardmedical.com](http://www.lombardmedical.com) clinician section.

## 10.0 Directions for Use

### 10.1. Pre-Procedure Planning

- Review CT angiograms (CTA) to assess aortic neck, aneurysm sac, and iliac artery anatomy.
- Select aortic and iliac endografts best suited for the anatomy.
- Confirm the recommended materials are available and ready for use.
- Anatomical requirements:
  - Iliofemoral access vessel size and morphology (minimal thrombus, calcium, and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.

- Proximal aortic neck lengths should be equal to or greater than of 15 mm with a lumen diameter of 18 mm to 28 mm, inclusive and a neck angle that does not exceed 60 degrees.
- Iliac artery distal fixation length should be equal to or greater than 15 mm with a lumen diameter of and 8 mm to 18 mm, inclusive.

**The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.**

### 10.2. General Use Information

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters, and guidewires should be employed during use of the Altura Endograft System.
- It is important to monitor blood loss throughout the procedure and if blood loss is excessive treat appropriately.

### 10.3. Patient Preparation

- Prepare and drape the patient for a surgical AAA procedure, in the event that the conversion to open repair is required.
- Position the patient on the imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
- Sedate the patient using standard practices.
- Anticoagulate the patient using standard practices. If patient is sensitive to heparin, an alternative anticoagulant should be considered.
- Access the left and right femoral arteries either percutaneously or surgically using the appropriate techniques.
- Establish adequate proximal and distal vascular control of both femoral vessels.

### 10.4. Implant Procedure: Aortic Endografts

#### 10.4.1. Implant Preparation

- Following standard techniques, introduce a guidewire and angiographic flush catheter into the aorta to assess the abdominal vasculature. Verify that no significant anatomical changes have occurred since the pre-operative CTA and the patient still meets the anatomical inclusion / exclusion criteria of the investigational plan. Verify the iliac and femoral arteries will allow endovascular access of introducer sheaths and catheters. Remove the angiographic flush catheter.
- Using standard technique, introduce bilateral 0.035" (0.89 mm) stiff guidewires. Advance both guidewires into the thoracic aorta.
- Introduce bilateral 14 F introducer sheaths.
- Carefully inspect the aortic delivery catheter packages for damage. Using sterile technique, carefully remove the aortic delivery catheters from the packaging and inspect the delivery catheters for damage. Do not use a delivery catheter if damage is present.
- Flush the delivery catheters with sterile heparinized saline through the guidewire flush port and implant flush port and while holding the tip elevated until a steady, constant stream exits the catheter.

**Note:** Saline may exit through the junction of the tip and top cap during flushing. If this occurs, the device may still be used.

- Rotate Aortic Delivery Catheter to ensure proper Top Cap rotation:
  - Hold the White Tip on one end and the Outer Sheath on the other end (away from the stent, **Figure 8**), and rotate the Top Cap by at least 180° (clockwise or counter clockwise, **Figure 9**) without pulling on the White Tip.
- **If the Top Cap does not rotate at least 180° with respect to the Outer Sheath, do not use the Aortic Delivery Catheter.**

#### 10.4.2. Tracking and Initial Positioning

**Note:** The following steps may be performed simultaneously with both aortic endografts or with one aortic endograft up to the point of final deployment. Ensure that a single aortic endograft is not fully implanted within the aorta without a paired aortic endograft positioned, aligned, and partially expanded alongside.

- Align the fluoroscopic guidance equipment to the preferred imaging orientation of the native anatomy.
- Load the aortic delivery catheters over the guidewires. Introduce the catheters and advance under fluoroscopic guidance so the middle septal radiopaque marker is positioned near the target renal artery. The catheters may be crossed. Care should be taken to maintain guidewire position during tracking.

**Note:** If the catheters are superimposed after being tracked into position, the delivery catheters may be advanced into the suprarenal aorta to position the catheters side by side. The image intensifier of the fluoroscope may also be rotated so that the two delivery catheters can be seen side by side.

- Under fluoroscopic guidance, rotate the delivery catheters so the septal radiopaque markers are parallel to and are inside of their respective guidewires. The endografts are rotated by turning the entire delivery catheter.
- Stabilize the front grip of the delivery catheters, ensuring that the catheters and the handles are straight.
- While stabilizing the front grip, slowly retract the outer sheath by rotating the rear grip of the delivery catheter handle as indicated on the handle. Deploy until the outer sheath radiopaque marker is aligned with the inner shaft radiopaque marker. Repeat this step for the contralateral device if simultaneously deploying both aortic endografts. To verify the position of the renal arteries, the outer shaft may be retracted so that the outer shaft radiopaque marker aligns with the first septal marker and then injected with diluted (50:50) contrast material through the implant flush port either manually or using a power injector. Do not exceed 300 psi while injecting. After the contrast injection, deploy until the outer sheath radiopaque marker is aligned with the inner shaft radiopaque marker.
- Confirm that the septal radiopaque markers are parallel to and are inside of their respective guidewires. The endografts are rotated by turning the entire delivery catheter.

**Note:** To move the stent expander to expand or collapse the endograft, the lock must be pressed. To secure the stent expander in place, release the lock.

- Stabilize the front grip and move the stent expander toward the rear of the delivery catheter to partially expand the exposed portion of the aortic endograft. Repeat this step for the contralateral device if simultaneously deploying both aortic endografts.

**Note:** If resistance is encountered during endograft rotation for septal alignment, one of the endografts may be advanced proximally or distally relative to the other to allow it to more freely rotate in the aorta.

- If necessary, rotate one or both aortic endografts so the flat surfaces appose each other. The endografts are rotated by turning the entire delivery catheter. If both aortic endografts are partially expanded and viewed so the endografts are side by side (generally an Anterior-Posterior orientation), rotate until the septal radiopaque markers are parallel to and are inside of their respective guidewires (**Figure 10**). If the endografts are not side by side, the delivery catheters may be rotated to properly align. To confirm correct alignment, the endografts may be viewed so they are superimposed on each other (generally a lateral orientation). In this view, the septal radiopaque markers of the two aortic endografts cross (**Figure 11**).

#### 10.4.3. Renal Positioning and Expansion

**Note:** The following steps may be performed simultaneously with both aortic endografts or with one aortic endograft up to the point of final deployment. Ensure that a single aortic endograft is not fully implanted within the aorta without a paired aortic endograft positioned, aligned, and partially expanded alongside.

- Verify the position of the renal arteries in relation to the renal radiopaque marker. Renal artery position can be determined by injecting diluted (50:50) contrast material through the implant flush port either manually or using a power injector. Do not exceed 300 psi while injecting.
- Carefully move one or both endografts so the graft material of each endograft is below the respective renal artery. Make certain graft material does not occlude the renal arteries. One or both aortic endografts may be collapsed to facilitate repositioning by advancing the stent expander forward.
- If one or both endografts were collapsed, re-expand the endografts to the partially expanded configuration using the stent expander. Confirm that the flat surfaces of the endografts are still properly apposed to each other. The flat surfaces of the endografts must appose each other before continuing with deployment.
- If simultaneously deploying both aortic endografts, stabilize the front grip of both delivery catheters and move the stent expanders toward the rear of the delivery catheter handle to expand the exposed portion of the aortic endografts. Expand until the endografts are well apposed to the vessel wall. If expanding only one aortic endograft at a time, expand the first endograft until the exposed portion of the endograft occupies a majority of the aortic lumen.
- Confirm correct renal positioning of the endografts using an injection of diluted (50:50) contrast material through the implant flush port. Do not exceed 300 psi while injecting. If repositioning is necessary, compress the exposed portion of one or both endografts by advancing the stent expander forward. After repositioning, re-expand the endografts and confirm the flat surfaces of the endografts appose each other.

**Note:** If one aortic endograft was being deployed, perform the tracking, positioning, alignment, and expansion steps of the contralateral aortic endograft. At this point of the procedure, proximal portions of both aortic endografts should be exposed, aligned, and expanded. Do not proceed in deployment until both endografts have reached this configuration.

- Remove the safety collar from the delivery catheter handle. Release the suprarenal stent by stabilizing the front grip and advancing the stent release forward until it snaps and locks together with the stent expander. Repeat for the contralateral aortic endograft.

**Note:** As the rear grip is rotated, the outer sheath is retracted and the inner member is advanced. This axially compresses and radially expands the endograft as it is unsheathed.

- This step may be done simultaneously with both aortic endografts or sequentially with each aortic endograft. While stabilizing the front grip, deploy the aortic endograft by rotating the rear grip of the delivery catheter handle as indicated on the handle. The aperture of the fluoroscope may be widened to allow for simultaneous visualization of the renal arteries and aortic bifurcation.

**Note:** Gentle compression on the delivery catheter handle may increase the expansion and flexibility of the aortic endograft. Use caution to avoid distal movement of the aortic endografts to avoid unplanned renal artery occlusion.

- Rotate the rear grips until the endografts are fully deployed.

#### 10.4.4. Delivery Catheter Removal

- Inspect the top cap region of the delivery catheters and ensure the radiopaque marker is positioned over the cylindrical portion of the leading collet.
- The delivery catheters should be removed one at a time. Stabilize the guidewire and introducer sheath positions. Under fluoroscopic guidance, carefully withdraw the delivery catheter over the guidewire. Take care not to catch the delivery catheter on the endograft or vascular introducer sheath as the catheter is removed. Gentle rotation of the delivery catheter may help prevent the catheter from catching on the endograft or vascular introducer sheath. Repeat this step with the contralateral delivery catheter.
- Leave the introducer sheaths and guidewires in place.

**Note:** The aortic endografts may be ballooned with compliant balloons prior to deployment of the iliac endografts. For instructions on ballooning, see section 10.6.

### 10.5. Implant Procedure: Iliac Endografts

#### 10.5.1. Implant Preparation

- Perform an angiogram of the iliac arteries to determine the location of the internal iliac arteries.

- Inspect the iliac delivery catheter packages for damage. Using sterile technique, carefully remove the iliac delivery catheters from the packaging and inspect the delivery catheters for damage. Do not use a delivery catheter if damage is present.
- Note:** The iliac delivery catheter has two small flush holes at the distal end of the top cap. The holes allow for the loaded endograft to be flushed. Heparinized saline should eject from these holes during flushing.
- Flush the implant and guidewire lumens of the delivery catheters with sterile heparinized saline while holding the tip elevated until a steady, constant stream exits the catheter. 10.5.2. Tracking and Initial Positioning
  - Align the fluoroscopic guidance equipment to the preferred imaging orientation of the native anatomy.
  - Under fluoroscopic guidance, advance the delivery catheter until the distal end of the iliac endograft is approximately 1 cm distal to the intended seal zone. Use care to ensure the iliac delivery catheter does not catch on the previously deployed aortic endograft. Stabilize the guidewire during tracking to prevent guidewire advancement.
  - Ensure the radiopaque marker of the top cap is outside the introducer sheath.

#### 10.5.3. Deployment

**Note:** The iliac delivery catheter deploys the iliac endograft by moving the top cap forward, releasing the distal portion of the endograft first.

- While stabilizing the front grip and ensuring that the catheters and the handles are straight, partially deploy the iliac endograft by rotating the rear grip of the delivery catheter handle as indicated on the handle.
- Internal iliac artery position can be determined by injecting diluted (50:50) contrast material through the implant flush port either manually or using a power injector. The top cap should be advanced so the top cap radiopaque marker is aligned with the distal edge of the endograft prior to injecting contrast in order to allow for adequate contrast flow. Do not exceed 300 psi while injecting.
- Minor repositioning of the iliac graft can be made proximally as the endograft makes initial contact with the iliac artery by advancing the delivery catheter forward.

**Note:** As the rear grip is rotated, the top cap is pushed forward and the leading collet is retracted. This axially compresses and radially expands the endograft as it is unsheathed. Gentle tension on the delivery catheter handle after the iliac endograft has made contact with the vasculature may increase the expansion and flexibility of the endograft. Do not reposition the iliac endograft distally once the endograft makes contact with the iliac artery.

- Release the iliac endograft by continuing to rotate the rear grip of the delivery catheter handle until the leading edge of the radiopaque marker of the top cap is aligned with the leading edge of the cylindrical portion of the collet.
- Confirm the proximal radiopaque markers on the iliac endograft have expanded off the leading collet.

#### 10.5.4. Delivery Catheter Removal

- Inspect the top cap region of the delivery catheter and ensure the radiopaque marker is aligned with the leading edge of the cylindrical portion of the leading collet.

**Note:** Gentle rotation of the delivery catheter may help prevent the catheter from catching on the endograft or vascular introducer sheath during removal. If resistance is felt during removal of the top cap through the introducer sheath, re-docking of the top cap may ease removal. To re-dock the top cap, rotate the rear grip in the opposite direction as indicated to return it to its pre-deployment position.

- Stabilize the guidewire and introducer sheath positions. Under fluoroscopic guidance, carefully withdraw the delivery catheter over the guidewire. Take care not to catch the delivery catheter on the endograft or vascular introducer sheath as the catheter is removed.
- Leave the introducer sheath and guidewire in place.
- If the iliac delivery catheter cannot be removed through the introducer sheath, remove the introducer sheath and delivery catheter together:
  - Prepare a new introducer sheath to be placed into the patient after the current introducer sheath is removed.
  - Retract the delivery catheter so that the proximal end of the iliac top cap is inside the distal segment of the iliac endograft. During retraction, gently rotate the delivery catheter to help prevent the catheter from catching on the endograft.
  - Rotate the rear grip of the iliac delivery catheter in the opposite direction as indicated on the handle. While rotating the rear grip, adjust the position of the delivery system so that the top cap continues to stay inside the iliac endograft with limited motion.
  - Continue to rotate rear grip until the rear grip can no longer be rotated. ○ Rotate rear grip ¼ turn in the direction of the arrows.
  - While maintaining guidewire position and hemostasis via manual pressure, remove the delivery catheter and introducer sheath together.
  - Insert the new introducer sheath into the patient over the indwelling guidewire.

10.5.5. If a minimum 20 mm overlap between the aortic and iliac endografts is not achieved, deploy an additional 13 mm iliac endograft over the junction of the previously deployed aortic and iliac endografts. Perform the iliac tracking, initial positioning, deployment, and removal steps for the additional iliac endograft.

10.5.6. Perform the iliac tracking, initial positioning, deployment, and removal steps on the contralateral side.

### 10.6. Endograft Balloon Dilatation

- Advance an appropriately sized compliant balloon catheter over the indwelling guidewires through both the left and right endografts. Advance both balloons to the aortic sealing zone of the aortic endografts.
- While using fluoroscopic guidance, carefully inflate both balloons simultaneously to appose the endografts to the aortic wall. Inflate the balloons by hand using partially filled 30 ml (cc) syringes.
- Deflate the balloons. If needed to fully open the suprarenal stent, advance the balloons proximally. Gently inflate both balloons simultaneously until the suprarenal stents are open.

- Deflate and reposition the balloons distally. Repeat the inflation procedure along the entire length of the endografts, including the iliac endografts.

#### **10.7. Post Implant Aortogram**

- Insert an angiographic flush catheter into the suprarenal aorta and perform an aortogram. Verify patency of renal and internal iliac arteries.
- Confirm there are no endoleaks or kinks and verify intended endograft placement.

#### **10.8. Procedure Completion**

- Close the vascular access sites using the appropriate techniques.
- Allow the patient to recover and be discharged accordingly.

#### **10.9. Disposal of delivery system device**

- After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy. If the device is damaged or the integrity of the sterilization barrier has been compromised, do not use the product and contact your Lombard Medical representative for return information.

### **11.0 Storage and Handling**

All components should be handled carefully to avoid damaging the device. Packages should be stored in a manner that protects the integrity of the package, the sterile barrier and prevents exposure to elevated temperatures. Store in a cool, dry place.

### **12.0 Warranty and Limitation of Liability**

Altura Medical® warrants that reasonable care was used in the design and manufacture of this product. Because Altura Medical has no control over the conditions of use, patient selection or handling of the device after it leaves its possession, Altura Medical does not warrant either a good effect or against an ill effect following its use. Altura Medical shall not be directly or indirectly responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. Altura Medical's sole responsibility, in the event Altura Medical determines the product was defective when shipped by Altura Medical, shall be the replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for use.

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Lombard Medical Ltd
Limited