

AORFLEX™

Delivery Device

The Flexibility To Deliver

Effectively Navigate Tortuous Anatomy


- Precise maneuvering within vessel walls

Advanced 1:1 Deployment Control

- Excellent torque, pushability and trackability
- High kink resistance and flexibility

The Confidence of Precise Positioning


- Platinum iridium marker delivers excellent radiopacity to accurately indicate position of the sheath during deployment



AORFLEX™

Delivery Device

Designed Specifically For Patients With Complex Anatomies



Federal (USA) law restricts this device to sale by or on the order of a physician

Refer to product IFU for full indications



15420 Laguna Canyon Road, Ste 260 | Irvine, CA 92618
Tel: 866 790 1447 | Fax: 949 379 3760 | www.lombardmedical.com
Email: USinfo@lombardmedical.com

© 2013 Lombard Medical Technologies Inc. The Lombard Medical logo, the Aorfix Endovascular Stent Graft logo, the Aorflex Delivery Device logo and the Aorfix and Aorflex names are trademarks or registered trademarks of Lombard Medical Technologies Inc.

PN002277