


Talos™ - Honeycomb Extension


Honeycomb Extension - Tapered 6mm

Product Code	Proximal Diameter (D1)	Distal Diameter (D2)	Length of covered graft (L1)	Length of Honeycomb section (L2)	Taper	Total Length (L1 + L2)	Profile (F)
TN2822-020140	28	22	20	140	6	160	20
TN3024-020140	30	24	20	140	6	160	20
TN3226-020140	32	26	20	140	6	160	20
TN3428-020140	34	28	20	140	6	160	20
TN3630-020140	36	30	20	140	6	160	20
TN3024-020180	30	24	20	180	6	200	20
TN3226-020180	32	26	20	180	6	200	20
TN3428-020180	34	28	20	180	6	200	20
TN3630-020180	36	30	20	180	6	200	20
TN4034-020180	40	34	20	180	6	200	22

Talos™ Thoracic Stent-Graft



Talos™ Honeycomb Extension



Lombard Medical Limited, Lombard Medical House, 4 Trident Park, Didcot, OX11 7HJ
 Tel: +44 (0)1235 750800 Email: globalinfo@lombardmedical.com Website: www.lombardmedical.com

Any information may be changed without prior notice. Talos™ is manufactured by MicroPort Endovastec and Distributed in Europe by Lombard Medical Limited.
 Not Available for sale in the United States. ©2026 Lombard Medical Limited. All rights reserved. D01504_ENG B



Talos™ -Thoracic Stent Graft

Fully Covered - Straight (mm)

Product Code	Proximal Diameter (D1)	Distal Diameter (D2)	Total Length (L1)	Profile (F)
T2828-160000	28	28	160	20
T3030-160000	30	30	160	20
T3232-160000	32	32	160	20
T3434-160000	34	34	160	20
T3636-160000	36	36	160	20
T3838-160000	38	38	160	22
T4040-160000	40	40	160	22
T4242-160000	42	42	160	22
T4444-160000	44	44	160	22
T2828-200000	28	28	200	20
T3030-200000	30	30	200	20
T3232-200000	32	32	200	20
T3434-200000	34	34	200	20
T3636-200000	36	36	200	20
T3838-200000	38	38	200	22
T4040-200000	40	40	200	22
T4242-200000	42	42	200	22
T4444-200000	44	44	200	22

Fully Covered - Tapered (mm)

Product Code	Proximal Diameter (D1)	Distal Diameter (D2)	Taper	Total Length (L1)	Profile (F)
T3026-160000	30	26	4	160	20
T3228-160000	32	28	4	160	20
T3430-160000	34	30	4	160	20
T3632-160000	36	32	4	160	20
T3834-160000	38	34	4	160	22
T4036-160000	40	36	4	160	22
T4238-160000	42	38	4	160	22
T4440-160000	44	40	4	160	22
T3026-200000	30	26	4	200	20
T3228-200000	32	28	4	200	20
T3430-200000	34	30	4	200	20
T3632-200000	36	32	4	200	20
T3834-200000	38	34	4	200	22
T4036-200000	40	36	4	200	22
T4238-200000	42	38	4	200	22
T4440-200000	44	40	4	200	22

Talos™ - Thoracic Stent Graft

Honeycomb Tapered (mm)

Product Code	Proximal Diameter (D1)	Distal Diameter (D2)	Length of covered graft (L1)	Length of honeycomb section (L2)	Taper	Total Length (L1 + L2)	Profile (F)
T2824-100060	28	24	100	60	4	160	20
T3026-100060	30	26	100	60	4	160	20
T3228-100060	32	28	100	60	4	160	20
T2620-160060	26	20	160	60	6	220	20
T2822-160060	28	22	160	60	6	220	20
T3024-160060	30	24	160	60	6	220	20
T3226-160060	32	26	160	60	6	220	20
T3428-160060	34	28	160	60	6	220	20
T3630-160060	36	30	160	60	6	220	20
T3832-160060	38	32	160	60	6	220	22
T4034-160060	40	34	160	60	6	220	22
T4236-160060	42	36	160	60	6	220	22
T4438-160060	44	38	160	60	6	220	22
T2824-160100	28	24	160	100	4	260	20
T3026-160100	30	26	160	100	4	260	20
T3228-160100	32	28	160	100	4	260	20
T3430-160100	34	30	160	100	4	260	20
T3632-160100	36	32	160	100	4	260	20
T3834-160100	38	34	160	100	4	260	22
T2620-200060	26	20	200	60	6	260	20
T2822-200060	28	22	200	60	6	260	20
T3024-200060	30	24	200	60	6	260	20
T3226-200060	32	26	200	60	6	260	20
T3428-200060	34	28	200	60	6	260	20
T3630-200060	36	30	200	60	6	260	20
T3832-200060	38	32	200	60	6	260	22
T4034-200060	40	34	200	60	6	260	22
T4236-200060	42	36	200	60	6	260	22

Talos™ - Honeycomb Extension

Honeycomb Extension - Straight (mm)

Product Code	Proximal Diameter (D1)	Distal Diameter (D2)	Length of covered graft (L1)	Length of Honeycomb section (L2)	Taper	Total Length (L1 + L2)	Profile (F)
TN2828-020100	28	28	20	100	0	120	20
TN3030-020100	30	30	20	100	0	120	20
TN3232-020100	32	32	20	100	0	120	20
TN3434-020100	34	34	20	100	0	120	20
TN3636-020100	36	36	20	100	0	120	20
TN2828-020140	28	28	20	140	0	160	20
TN3030-020140	30	30	20	140	0	160	20
TN3232-020140	32	32	20	140	0	160	20
TN3434-020140	34	34	20	140	0	160	20
TN3636-020140	36	36	20	140	0	160	20
TN2424-020200	24	24	20	200	0	220	20
TN2626-020200	26	26	20	200	0	220	20
TN2828-020200	28	28	20	200	0	220	20
TN3030-020200	30	30	20	200	0	220	20

Honeycomb Extension - Tapered 4mm

Product Code	Proximal Diameter (D1)	Distal Diameter (D2)	Length of covered graft (L1)	Length of Honeycomb section (L2)	Taper	Total Length (L1 + L2)	Profile (F)
TN2824-020100	28	24	20	100	4	120	20
TN3026-020100	30	26	20	100	4	120	20
TN3228-020100	32	28	20	100	4	120	20
TN3430-020100	34	30	20	100	4	120	20
TN3632-020100	36	32	20	100	4	120	20
TN4036-020100	40	36	20	100	4	120	22
TN2824-020140	28	24	20	140	4	160	20
TN3026-020140	30	26	20	140	4	160	20
TN3228-020140	32	28	20	140	4	160	20
TN3430-020140	34	30	20	140	4	160	20
TN3632-020140	36	32	20	140	4	160	20