



AlturaTM

Abdominal Aortic Stent-Graft System


Revolutionary Design, Extraordinary Advantages

Simplified Delivery, Precise Placement




Unique Bilateral D-Stent Design:

- Braided nitinol stent provides exceptional flexibility and vessel conformability
- Two sets of suprarenal anchors on each aortic stent
- No contralateral limb cannulation
- Radiopaque markers for precise positioning and easy septum alignment



Ultra-Low Profile 14F Delivery System:

- Designed to improve access and minimise vessel trauma
- Intuitive five-step delivery system
- Streamlined endograft delivery ensures controlled and accurate placement
- Built-in contrast injection ports simplify the procedure




Maximises Sealing Zone:

- Targets each vessel and provides four points of accuracy
- Ability to independently position and offset aortic endografts maximises seal in infrarenal neck

Traditional



Altura™




Repositionable for Accurate Placement:

- Aortic endografts may be repositioned for precise infrarenal placement

Retrograde Iliac Deployment:

- Unique one-step retrograde iliac endograft deployment, for accurate placement at the iliac landing zone




“Altura’s unique design eliminates gate cannulation, which can be challenging and time-consuming with a traditional stent. Coupled with the retrograde iliac deployment, these innovations make my procedures predictable, extremely precise, and easier.”


Prof. Dainis Krievins
Vascular Surgeon, Director, Institute of Science, Stradins University Hospital, Riga, Latvia

Product and Ordering Information

Aortic Endograft (mm)



Iliac Endograft (mm)



The Altura™ Endograft portfolio consists of six product codes, which allows the treatment of a broad range of patients, with simplified case planning and inventory management.

Aortic Endograft (mm)

Aortic Treatment Range*	Endograft Diameter (D1)	Suprarenal Stent Diameter (D2)	Product Code
18-22	24	29	01-AA-24-090
21-25	27	33	01-AA-27-090
24-28	30	36	01-AA-30-090

Iliac Endograft (mm)

Iliac Treatment Range*	Distal Diameter (D3)	Product Code
8-11	13	01-AI-13-065
11-15	17	01-AI-17-065
15-18	21	01-AI-21-065

*Treatment range measured from inner lumen



Lombard Medical Limited, Lombard Medical House, 4 Trident Park, Didcot, OX11 7HJ
Tel: +44 (0)1235 750800 Email: globalinfo@endovastec-lm.com Website: www.lombardmedical.com

All cited trademarks are the property of their respective owners in respective countries. These devices should only be used by a licenced physician. Indicators, contraindications, warnings, potential adverse events and directions for use can be found in Instructions for Use and product labelling supplied with each device. The promotion, distribution or sale of these devices is prohibited in countries where the registration certificates or any other licences have not been obtained. This brochure is for the use of the individuals to whom it has been provided and is not to be distributed widely as an advertisement. Any information may be changed without prior notice. Not Available for sale in the United States.

©2026 Lombard Medical Limited. All rights reserved. D01513_ENG A