




*Minos*<sup>TM</sup>

Abdominal Aortic  
Stent-Graft and Delivery System

Expanding the Possibilities of EVAR with  
Ultra-Low Profile Design

Reliable Design, Safety and Durability Ensured



**Tri-Modular Design:**  
- Adjustable overlapping length  
- Reduced inventory

**Hydrophilic Coating:**  
- Supports smooth introduction

**Ultra-Low Profile:**  
- 12F, 14F and 16F




**Helical Limb Stent:**  
- Highly flexible  
- Reduces the risk of occlusion

**Dual Proximal Stent:**  
- Enhances conformability and proximal seal

**Seamless PET Graft:**  
- Durable multi-filament weave

**Laser Cut Bare Stent Hooks:**  
- Provide active fixation

**Radiopaque Markers**



12 distinctive markers support accurate positioning and excellent visualisation



# Ultra-Low Profile Design, Expanding the Possibilities for Treating Challenging Anatomies

## Product and Ordering Information

### Main Body Stent-Graft


The Main Body Stent-Graft	CM22-100	CM24-100	CM26-100	CM28-100	CM30-100	CM32-100	CM34-100
Outer Sheath Profile (OD)	14F			16F			
Delivery System Working Length	550mm						

Example: CM22-100

CM = Main Body

D1 = 22mm Main Body proximal diameter

L1 = 100mm Ipsilateral graft length



### Iliac Stent-Graft



Iliac Stent-Graft	CL10-80	CL13-80	CL16-80	CL18-80	CL20-80	CL24-80
	CL10-100	CL13-100	CL16-100	CL18-100	CL20-100	CL24-100
	CL10-120	CL13-120	CL16-120	CL18-120	CL20-120	CL24-120
	CL10-140	CL13-140	CL16-140	CL18-140	CL20-140	CL24-140
Outer Sheath Profile (OD)	12F					
Delivery System Working Length	610mm					

Example: CL13-80

CL = Iliac Stent-Graft

D3 = 13mm Distal Limb diameter

L3 = 80mm Iliac Stent-Graft Length



**Lombard Medical Limited, Lombard Medical House, 4 Trident Park, Didcot, OX11 7HJ**  
**Tel: +44 (0)1235 750800 Email: [globalinfo@lombardmedical.com](mailto:globalinfo@lombardmedical.com) Website: [www.lombardmedical.com](http://www.lombardmedical.com)**

All cited trademarks are the property of their respective owners in respective countries.

These devices should only be used when prescribed by a licenced physician. Indicators, contraindications, warnings, potential adverse events and directions for use can be found in Instructions for Use and product labelling supplied with each device.

The promotion, distribution or sale of these devices is prohibited in countries where the registration certificates or any other licences have not been obtained. This brochure is not intended for use in aforesaid countries. This brochure is only delivered to specific people for their reference, is not to be distributed widely as an advertisement. Any information may be changed without prior notice.

**Not Available for sale in the United States.**

©2020 Lombard Medical Limited. All rights reserved. MP00004. Version 2.0